

Manual

EconoTaq PLUS 2X Master Mix

For Research Use Only. Not for use in diagnostic procedures.

IMPORTANT
-20 °C storage required
immediately upon receipt

Manual

EconoTaq PLUS 2X Master Mix

Contents

1. Product description	3
2. Introduction to PCR	3
3. Product designations and kit components	4
4. Storage conditions	4
5. PCR setup	4
6. PCR cycling conditions	5
7. PCR guidelines	5
8. Technical support and product guarantee	6

Manual

EconoTaq PLUS 2X Master Mix

1. Product description

EconoTaq™ PLUS 2X Master Mix is a ready-to-use PCR master mix containing agarose gel loading buffer. A PCR reaction is set up simply by combining the EconoTaq PLUS Master Mix with template DNA, primers and water. After the reaction is complete, it can be loaded directly onto an agarose gel. All other necessary components are provided in the Master Mix, which contains the following: 0.1 U/μL of EconoTaq DNA Polymerase, Reaction Buffer (pH 9.0), 400 μM dATP, 400 μM dGTP, 400 μM dCTP, 400 μM dTTP, 3 mM MgCl₂ and a proprietary PCR Enhancer/Stabilizer. For difficult templates, such as those containing GC content up to 75%, it has been found that EconoTaq PLUS Master Mix performs well with denaturation temperatures up to 98 °C.

PCR activity: EconoTaq 2X Master mix is tested in DNA amplification using a variety of templates and primers.

Activity determination: One unit of EconoTaq DNA Polymerase catalyses the incorporation of 10 nmoles of dNTP into acid-insoluble material in 30 minutes at 70 °C in 50 mM Tris-HCl (pH 9.0), 50 mM NaCl, 5 mM MgCl₂, 200 μM dGTP, dATP, dTTP, dCTP (a mix of unlabeled and [³³P]dCTP), 10 μg Activated Calf Thymus DNA and 0.1 mg/mL BSA.

Absence of endonuclease or nicking activity: Incubation of 10 U of EconoTaq DNA Polymerase with 1 μg of supercoiled pBR322 DNA for 16 hours at 70 °C results in no detectable conversion to relaxed or linear forms detectable by agarose gel electrophoresis.

Absence of exonuclease activity: Incubation of 10 U of EconoTaq DNA Polymerase with 1 μg of *Hind*III-cut lambda DNA for 16 hours at 70 °C resulted in no smearing of bands on agarose gels.

Purity: EconoTaq DNA Polymerase is >99% pure as determined by SDS-PAGE. There is no detectable DNA contamination.

2. Introduction to PCR

The Polymerase Chain Reaction (PCR) is a powerful technique that amplifies specific DNA sequences using multiple cycles of a 3-step process. The first step involves denaturation of a double-stranded DNA template at a high temperature. In the second step, sequence-specific primers anneal to complementary sites flanking the target sequence. In the third step, a thermostable DNA polymerase extends the annealed primers, thereby copying the original DNA sequence. The newly synthesised DNA becomes the template for subsequent DNA amplification, doubling the amount of template with each cycle. These steps are typically repeated 25-35 times, resulting in a 10⁵-10⁹ fold increase in the amount of target DNA.

Manual

EconoTaq PLUS 2X Master Mix

3. Product designations and kit components


Product	Kit size	Catalogue number	Reagent description	Part number	Volume
EconoTaq PLUS 2X Master Mix*	50 rxns**	30035-0	EconoTaq PLUS 2X Master Mix	F93482-0	1.25 mL
	500 rxns	30035-1	EconoTaq PLUS 2X Master Mix	F93482-1	12.5 mL
	1,000 rxns	30035-2	EconoTaq PLUS 2X Master Mix	F93482-1	25.0 mL

* EconoTaq PLUS 2X Master Mix contains all of the components necessary to perform PCR amplification, except for template and primers. The Master Mix also contains agarose gel loading buffer..

** Trial size.

4. Storage conditions

EconoTaq PLUS 2X Master Mix can be stored at -20 °C for 12 months. For everyday use, an aliquot can be stored at +4 °C for up to three months. The EconoTaq PLUS 2X Master Mix needs to be mixed well prior to use. It is stable for ten freeze-thaw cycles.



5. PCR setup

5.1. **Materials supplied by the user.** PCR amplification is performed by adding template DNA, primers and water to the EconoTaq PLUS 2X Master Mix. The following components must be supplied by the user:

- Template DNA (10-50 ng of plasmid DNA; 50-200 ng of genomic DNA)
- Forward primer (100 µM)
- Reverse primer (100 µM)
- Nuclease-free water
- Thermal cycling apparatus

Manual

EconoTaq PLUS 2X Master Mix

5.2. **Reaction setup.** Set up PCR amplifications of the desired size, according to the following:

				Final concentration
EconoTaq PLUS 2X Master Mix	12.5 µL	25.0 µL	50.0 µL	1 X
Forward primer (100 µM)	0.25 µL	0.5 µL	1.0 µL	1 µM
Reverse primer (100 µM)	0.25 µL	0.5 µL	1.0 µL	1 µM
DNA template (10 ng/µL)	1.0 µL	1.0 µL	1.0 µL	
Water, nuclease-free	11.0 µL	23.0 µL	47.0 µL	
Total reaction volume	25.0 µL	50.0 µL	100.0 µL	

5.2. Gently mix the PCR components in a thin-walled reaction tube and spin briefly in a microcentrifuge. Add a drop of mineral oil if the thermal cycler does not have a heated lid.

6. PCR cycling conditions

6.1. Pre-heat the thermal cycler to 94 °C.

6.2. For initial denaturation of target template DNA, incubate the reactions at 94 °C for two minutes.

6.3. Denature, anneal and extend the DNA according to the following for subsequent cycles of amplification:

Cycling step	Temperature	Time	# of Cycles
Initial denaturation*	94 °C	2 min	1
Denaturation*	94 °C	15-30 sec	25-35
Annealing**	50-65 °C	15-30 sec	
Extension	72 °C	1 min/kb	
Final extension	72 °C	5-10 min	1
Hold	4 °C	Indefinitely	1

* EconoTaq PLUS performs well with denaturation temperatures up to 98 °C.

** Anneal at T_m of primer \pm 2 °C. See PCR guidelines section.

6.4. After completion of the PCR, a 5 µL aliquot of the reaction can be loaded directly onto an agarose gel for analysis or size selection.

Manual

EconoTaq PLUS 2X Master Mix

7. PCR guidelines

Optimisation of PCR conditions is required for amplification involving templates with high GC content, internal secondary structure or products greater than 5 kb. The following guidelines can be used to improve the success of amplifying these templates:

- 1) **Cold reaction setup.** Reactions using EconoTaq PLUS Master Mix yield best results when set up on ice and maintained at 4 °C prior to amplification. The EconoTaq Polymerase has residual activity at temperatures above 4 °C that can cause non-specific background amplification. Add primers last to reactions just prior to target addition and incubation/cycling.
- 2) **Template DNA.** Use of purified, high quality template DNA enhances the success rate of PCR. We recommend using 10-50 ng of plasmid DNA and 50-200 ng of genomic DNA. Template should be suspended in water, rather than EDTA-containing solution such as TE buffer.
- 3) **Primer design.** Oligonucleotide primers for PCR should be 20-25 bases in length and should have a GC content of 40 - 60%, with the GC bases evenly spaced in the primer. Self-annealing of primers leads to production of primer-dimers, which can diminish the amount of authentic product; therefore, the 3' end of each primer should not be complementary to itself or to the 3' end of the opposing primer. The melting temperature (T_m) of the primers should be within 5 °C of each other. The final concentration of each primer in the reaction should be 1 μ M.
- 4) **Annealing.** PCR enhancers in EconoTaq PLUS Master Mix may lower primer T_m slightly. Optimise the reaction conditions by performing the reaction starting with an annealing temperature 2 °C below the calculated melting temperature of the primers. Vary the annealing temperature to optimise for the primer/template combination. Reducing the annealing temperature improves yield and reaction efficiency, while increasing annealing temperature improves specificity.
- 5) **Denaturation temperature.** EconoTaq PLUS 2X Master Mix will amplify most routine DNA templates using a denaturation temperature of 94 °C. For GC-rich templates, with GC content up to 75%, a denaturation temperature of 98 °C has been found to improve results.

8. Technical support and product guarantee

If you require any further support, please do not hesitate to contact our Technical Support Team: techsupport@lgcgroup.com.

Product guarantee: LGC Biosearch Technologies™ guarantees that this product will perform as specified for one year from the date of shipment. Please avoid using reagents for greater than one year from receipt.



   @LGCBiosearch | biosearchtech.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2023. All rights reserved. GEN/953/EK/0123

For Research Use Only. Not for use in diagnostic procedures.



BIOSEARCH™
TECHNOLOGIES

GENOMIC ANALYSIS BY LGC