

# BigEasy® v2.0 Linear Cloning System

# Maximum insert stability

- · Clone "unclonable" DNA with this unique, high-efficiency system
- Stabilize up to 30 kb of difficult DNA sequence
- Create bias-free libraries from A/T-rich or G/C-rich genomes
- · Clone gene clusters or operons

BigEasy Systems containing the unique pJAZZ vector are ideal for cloning large DNA fragments, and for your most difficult DNA sequences that do not clone in other cloning vectors. Choose BigEasy for cloning single inserts up to 30 kb, or constructing bias-free, large-insert genomic libraries. The novel design of the pJAZZ vector includes multiple functional advantages over other vectors.

- pJAZZ is maintained as a linear molecule, so the vector does not supercoil. The ends of the vector can rotate freely (Fig. 1) reducing torsional stress.
- pJAZZ is maintained at low copy number (5-10/cell) in BigEasy-TSA™ Electrocompetent Cells, which are required for transformation and propagation. The copy number can be induced 5-10X for efficient DNA recovery and the vector can be isolated with standard plasmid prep methods.
- pJAZZ incorporates Lucigen's CloneSmart® technology for transcription-free cloning, which further increases insert stabilit.y

These unique advantages allow BigEasy Systems to stabilize large, repetitive, A/T-rich and G/C-rich or otherwise unclonable sequences (Fig. 2 and 3).





Figure 1. pJAZZ-OC and pJAZZ-OK linear vectors. RepA, replication factor and low copy origin of replication (~5-10 per cell; inducible 5-10 fold); Camr - chloramphenicol resistance gene; Kanr - kanamycin resistant gene; telN - protelomerase gene; cB - replication regulator. Approximate positions of transcription terminators (T) are indicated.



vector. This DNA was unclonable in all other vectors.

#### 96% AT rich DNA




Figure 2b. Sequence trace of a Piromyces clone showing extremely high AT content (96%).

#### 6-20 kb Oxytricha inserts




Figure 3. Oxytricha trifallax genomic DNA (75-85% AT) was sheared to 6-20 kb and cloned into the pJAZZ linear vector. NotI digests of mini-prep DNA are shown

**Lucigen Corporation** 608 831 9011

Toll Free: 888 575 9695 608 831 9012

www.lucigen.com custserv@lucigen.com techsupport@lucigen.com





## **CLONING KITS & VECTORS**

## lucigen.com/cloning



## BigEasy® v2.0 Linear Cloning System cont.

Products	Size	Cat. No.	Price
BigEasy v2.0 Linear Cloning Kit (pJAZZ-OC Blunt Vector) w/BigEasy-TSA Electrocompetent Cells (SOLOs)	5 reactions	43018-1	\$330
	10 reactions	43018-2	\$564
	20 reactions	43018-3	\$1010
BigEasy v2.0 Linear Cloning Kit (pJAZZ-OC Notl Vector) w/BigEasy-TSA Electrocompetent Cells (SOLOs)	5 reactions	43024-1	\$330
	10 reactions	43024-2	\$564
	20 reactions	43024-3	\$1,010
BigEasy v2.0 Linear Cloning Kit (pJAZZ-OK Blunt Vector) w/BigEasy-TSA Electrocompetent Cells (SOLOs)	5 reactions	43036-1	\$330
	10 reactions	43036-2	\$564
	20 reactions	43036-3	\$1010
BigEasy v2.0 Linear Cloning Kit (pJAZZ-OK Notl Vector) w/BigEasy-TSA Electrocompetent Cells (SOLOs)	5 reactions	43042-1	\$330
	10 reactions	43042-2	\$564
	20 reactions	43042-3	\$1010
BigEasy-TSA Electrocompetent Cells	6 reactions	60224-1	\$164
(≥ 4 × 10 <sup>10</sup> cfu/μg) (SOLOs)	12 reactions	60224-2	\$278
	24 reactions	60224-3	\$514

### ORDER INFORMATION

The BigEasy® Linear Cloning Kit includes: Dephosphorylated pJAZZ® Vector pre-cut at either a *Smal* (blunt) or *Not*1 site, CloneSmart® DNA Ligase, CloneDirect® 10X Ligation Buffer (includes ATP), DNATerminator® End Repair Enzyme & 5X End Repair Buffer (Blunt Kit only), Sequencing Primers, Positive Control Insert DNA, BigEasy-TSA Electrocompetent Cells in SOLO packaging (1 transformation per tube), Recovery Medium, Transformation Control DNA, and complete protocols. BigEasy-TSA Electrocompetent Cells are also available separately.

 Lucigen Corporation

 Phone:
 608 831 9011

 Toll Free:
 888 575 9695

 Fax:
 608 831 9012

bsi. ISO
13485
Medical Devices
Quality
Management

