

  
**Lucigen**<sup>®</sup>

Simplifying Genomics

epicentre<sup>®</sup>

## RNase I

Cat. No. N6901K

## 1. Introduction

RNase I preferentially degrades single-stranded RNA to individual nucleoside 3' monophosphates by cleaving every phosphodiester bond.<sup>1</sup> By comparison, other ribonucleases cleave only after specific residues (e.g., RNase A cleaves 3' to pyrimidine residues). Thus, RNase I is useful for removing RNA from DNA preparations,<sup>2</sup> detecting mismatches in RNA:RNA and RNA:DNA hybrids<sup>2,3</sup>, and analyzing and quantifying RNA in ribonuclease protection assays (RPA).<sup>4,5</sup> The enzyme is completely inactivated by heating at 70°C for 20 minutes in the presence of 5 mM dithiothreitol (DTT), eliminating the requirement to remove the enzyme prior to many subsequent procedures.

## 2. Product Designations and Kit Components

Product	Kit Size	Catalog Number	Reagent Description	Part Numbers	Volume
RNase I, <i>E. coli</i>	1,000 Units	N6901K	RNase I (10 U/μL)	E0067-10D1	100 μL
			DTT (0.1 M)	SS000065-D1	2.5 mL
			RNase I Dilution Buffer	SS000255-D1	1 mL
			10X TNE Buffer	SS000806-D1	5 mL

## 3. Product Specifications

**Storage:** Store only at -20°C in a freezer without a defrost cycle.

**Storage Buffer:** RNase I is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl, and 0.1 mM EDTA.

**RNase I Dilution Buffer:** A 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 0.1 M NaCl, and 0.1 mM EDTA.

**10X TNE Buffer:** 100 mM Tris-HCl (pH 7.5), 1 M NaCl, and 10 mM EDTA.

**Unit Definition:** One unit degrades 100 ng of *E. coli* ribosomal RNA per second into acid-soluble nucleotides at 37°C.

**Quality Control:** RNase I is function-tested in a reaction containing 10 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1 mM EDTA, and 60 μg of *E. coli* ribosomal RNA with varying amounts of enzyme.<sup>6</sup>

**Contaminating Activity Assays:** RNase I is free of detectable exo- and endodeoxyribonuclease activities as judged by incubation of 1 μg of various DNA substrates with 4 x 10<sup>6</sup> U of enzyme at 37°C for 16 hours.

## 4. Protocol for Removing RNA from DNA Preparations

RNase I can be used in place of RNase A for removing RNA from DNA preparations. In contrast to RNase A, RNase I effectively degrades contaminating RNA to mono- and dinucleotides that will not interfere with visualization of small DNA molecules. After RNA removal, the enzyme can be inactivated by heating at 70°C for 20 minutes in the presence of 5 mM DTT.

## Protocol

1. Isolate DNA from 1-2 mL of overnight bacterial culture using a standard alkaline lysis procedure.<sup>5</sup>
2. After ethanol precipitation, suspend the DNA in 1X TNE buffer (page 2) at a concentration appropriate for subsequent applications (see Notes II below).
3. Dilute RNase I ten-fold with RNase I Dilution Buffer and add 1.5-2 U to the DNA preparation.
4. Incubate at 37°C for 30 minutes to degrade contaminating RNA.
5. Add DTT to a final concentration of 5-10 mM.
6. Incubate at 70°C for 20 minutes to inactivate the enzyme.

## Notes

**Reaction Buffer:** Incubation with RNase I can be performed simultaneously with the digestion of plasmid DNA by restriction endonucleases. RNase I maintains  $\geq 90\%$  activity in buffers containing between 100 mM to 200 mM salt (either NaCl or KOAc). The activity of the enzyme is also relatively constant over a pH range of 7.0-8.8. Therefore, if the restriction endonuclease buffer is within these parameters, RNase I digestion can be performed in the restriction endonuclease buffer.

**Enzyme Dilution:** Diluted enzyme may be stored for up to two months at  $-20^{\circ}\text{C}$  in a freezer without a defrost cycle.

## 5. References

1. Shen, V. and Schlessinger, D. (1982) *The Enzymes XV (Part B)*, 501.
2. Winter, E. *et al.*, (1985) *Proc. Natl. Acad. Sci. USA* **82**, 7575.
3. Myers, R.M. *et al.*, (1985) *Science* **230**, 1242.
4. Sambrook, J. *et al.*, (1989) in: *Molecular Cloning: A Laboratory Manual (2nd ed.)*, Cold Spring Harbor Laboratory Press, New York.
5. Saccomanno, C.F. *et al.*, (1992) *BioTechniques* **13**, 847.
6. Corbishley, T.P. *et al.*, (1984) *Meth. Enzymatic Anal.* **4**, 134.

*Epicentre is a trademark of Illumina, Inc. and/or its affiliate(s) in the U.S. and other countries, and is used under license.*