


Lucigen[®]

Simplifying Genomics

epicentre[®]

Nuclease-Free Water

Cat. No. W7350ML

Nuclease-Free Water is useful for many molecular biological and biochemical applications.

1. Product Designations and Kit Components

| Product | Kit Size | Catalog Number | Reagent Description | Part Numbers | Volume |
|---------------------|----------|----------------|---------------------|--------------|--------|
| Nuclease-Free Water | 50 mL | W7350ML | Nuclease-Free Water | SS000772-D1 | 50 mL |

2. Product Specifications

Storage: Store at room temperature.

Molecular Weight: 18.016 daltons.

Contaminating Activity Assays: Nuclease-Free Water is prepared by autoclaving deionized water. The water is free of detectable exo- and endonuclease and RNase activities.

Epicentre is a trademark of Illumina, Inc. and/or its affiliate(s) in the U.S. and other countries, and is used under license.